




## TruSteel™

### TruSteel Infusion set

When using the TruSteel infusion sets, there is no need to “Fill Cannula” in the load sequence.

Your skin must be pinched into a skin fold before inserting the TruSteel Infusion Set.



### MONTHLY CONSUMABLE ORDERS MADE EASY

Request a recurring standard monthly delivery of your Tandem and Dexcom consumables from the following pharmacies:



012 426 4535  
diabetes@medipost.co.za



011 712 6000  
063 101 0911 (WhatsApp)  
brianm@cdediabetes.co.za

Our Customer Service Team remain available to you via the Customer Service number 0861 339 266 for any Tandem insulin pump related concerns or troubleshooting

Remember to provide the product name, quantity and NAPPI code to the pharmacy.