



t:slim X2™
Insulin Pump

How much insulin do you need to put into your t:slim cartridge?

When preparing your t:slim insulin cartridge there is a minimum of **120 units** of insulin required. and a maximum of **300 units**.

On entering the “Fill Tubing” step in the load process, there is a minimum tubing fill of **10.2 units** of insulin.

Once you have observed 3 droplets of insulin at the end of your infusion set, you may stop the fill tubing process and continue with the next step.



[Click here to see how to change a cartridge](#)

Our Customer Service Team remain available to you via the Customer Service number 0861 339 266 for any Tandem insulin pump related concerns or troubleshooting

MONTHLY CONSUMABLE ORDERS MADE EASY

Request a recurring standard monthly delivery of your Tandem and Dexcom consumables from the following pharmacies:

Medipost
Pharmacy

012 426 4535
diabetes@medipost.co.za

CDE Cardiovascular
Diabetes
Education
Your value driven partner in care

011 712 6000
063 101 0911 (WhatsApp)
brianm@cdediabetes.co.za

Remember to provide the product name, quantity and NAPPI code to the pharmacy.